

Project for Course Φ468: Light and Matter interactions and optical constants of Solids

Student ID Number: ph5207

Student name: Theodoros Laoumzis

Terahertz Time Domain Spectroscopy parameter analysis

The goal of this project was to extract the complex refractive index of a 2mm material. A set of measurements was given. The first data file is containing the reference and the second containing the sample measurements. Using programming language Python, I tried to write a code for the data analysis of the project.

Process of data analysis

First, I uploaded the data files and then plotted the two pulses. Then, I passed the measurements to frequency, performing a Fourier Transform (FFT) and plotted them again to see the results.

Afterwards, I divided the two signals to get the experimental Transfer Function ($T(f)$). I plotted the transfer function's real and imaginary parts to see the relation between them, then plotted the natural logarithm of $T(f)$ versus frequency.

To fit the theoretical transfer function, we assume that $n \gg k$, so that the theoretical (zero reflection) $T(f)_{th}$ is:

$$H_{zero}(\omega) = \frac{T_{10}(\omega) \cdot P_1(\omega, d) \cdot T_{01}(\omega)}{P_0(\omega, d)}$$

where T_{10} and T_{01} are the transmission coefficients.

Unwrapping the phase of the transfer function allowed me to calculate the real part of the refractive index: $n(\omega) = 1 + \frac{\varphi(\omega)c}{\omega d}$, where

- i. $\varphi(\omega)$ is the phase of the transfer function
- ii. $c = 3 \cdot 10^8 \text{ m/s}$ is the speed of light
- iii. ω is the frequency and
- iv. $d = 0.002 \text{ m}$ is the sample material's thickness.

I also calculated and plotted the absorbance $a(\omega)$ and extinction coefficient, k .

Finally, using the Newton – Raphson method, I tried to fit the theoretical transfer function to the experimental one and extract the complex refractive index, N.

The equations I used for the fitting process are the following:

- $F(\omega) = \ln\left(\frac{4\bar{n}_0\bar{n}_1}{(\bar{n}_0 + \bar{n}_1)^2}\right) - \frac{i\omega d(\bar{n}_1 - \bar{n}_0)}{c}$
- $dF = \frac{1}{\bar{n}_1} - \frac{2}{\bar{n}_0 + \bar{n}_1} - \frac{i\omega d}{c},$

where \bar{n}_0 and \bar{n}_1 are the complex refractive indices of the air and sample material.